Unravelling the little cherry disease complex at European scale to improve transnational diagnostics and management of the disease (EURAVELCH)



K. De Jonghe

Flanders Research Institute for Agricultural, Fisheries and Food (ILVO)



Aim of the project

Getting a better insight in the pathogenesis of the host plants with little cherry disease symptoms



Project partners

AGES: Thomas Leichtfried AGROSCOPE: Jean Sébastien Reynard ANSES: Jean Philippe Renvoisé CRA-W: Thibaut Olivier CTIFL: Julie Castaing ILVO: Kris De Jonghe JKI: Wilhelm Jelkmann NAKTUINBOUW: Ellis Meekes ULg: Sébastien Massart

by means of high throughput sequencing (NGS) (unravelling the little cherry disease complex)

- Main project tasks

- 1. Obtaining more information about the EU-wide spread of little cherry viruses (LChV), and raising awareness.
- 2. Optimization and standardization of diagnostics, taking into account the population diversity at national, European and worldwide levels (Networking & NGS Method performance study).
- 3. Enhancing knowledge on the little cherry disease epidemiology through NGS.
- 4. Reflecting on consequences for future EU quarantine legislation.

Case study: Epidemiology/diagnostics on a Prunus avium *rootstock carrying 3 ≠ cv. grafts*

Link disease symptoms and



targeted diagnostics:

- Several different specific PCR methods
- ✓ PCR result interpretation depends on
 - Choice of tests
 - Seasonal variation
 - Reliability of the individual tests
- ✓ Sequencing results:
 - Fragmented
 - Low resolution

A : cv. Hedelfinger
B : cv. Bigarreau-Burlat
C : cv. Bigarreau-Reverchon
Grafted 5-6y ago on same rootstock

rootstock (wild P. avium) CVA, PDV, LChV1, ? (not sequenced by NGS yet) untargeted diagnostics (NGS):

- Unbiased and broader view
- $\checkmark\,$ NGS result interpretation depends on
 - Sequencing strategy

 Virus enrichment
 - dsRNA, smallRNA, totRNA
 - Depth of sequencing
 - Bioinformatics pipeline
- ✓ Sequencing results:
 - Across genome variation study possible

* Admera Health, USA

2016-A-198

<u>Case study strategy</u>: SmallRNA sequencing NEBNext Ultra RNA Lib kit * Illumina NextSeq v2 * 3M reads; 1x 50bp VirusDetect automated pipeline SearchSmallRNA ref. genome mapping

Direct advantages of the NGS strategy:

- > Strawberry latent ringspot virus was not picked up by traditional strategy !
- Detection of novel viruses.
- Endpoint PCR detection of LChV-1&2 often erratic.
- Full genomes of CVA, PDV, LChV1, LChV2 obtained.

Allows standardisation of diagnostic methods

Better epidemiological insight



A: CVA, PDV + novel virus ?
B: CVA, PDV, SLRSV, LChV1, LChV2
C: CVA, PDV, SLRSV, LChV1

? ≠ virus/host interaction ? (cv ? rootstock ? ...)

Acknowledgements: Yoika Foucart (lab & data analysis) Annelies Haegeman (bioinformatics support)





RI 16/A-198 - EURAVELCH federal public service HEALTH, FOOD CHAIN SAFETY AND ENVIRONMENT

